IT26-1YS

Trazanto

Indicador Integrador de Un Punto

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Trazanto

Unique point Integrator Indicator

For Steam sterilization processes.

Usage

Control of Steam sterilization processes between 121 °C and 135 °C.

Applicable Regulation

Designed under Quality Management Systems standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11140-1:2014.

Classification

Class 1, according to risk.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Characteristics

Type 5 Integrator Indicator.

One-side-laminated strips printed with indicator ink.

Size: 70 mm x 23 mm (approximate).

Initial color: **yellow**.

Final color: **black**.

100 % Toxic Heavy Metals free

Integration condition calibrated with the kill time of a 10^6 *Geobacillus stearothermophilus* spores population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

The product ensures an adequate control of the efficacy of sterilization processes (temperature, time, quality of steam).

To interpret results visually, use the reference found in the indicator. Optionally, the readout of the chemical indicator can be done using Trazanto Automatic Quality Control and Traceability System (see instructions for use in device's manual).

Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

Storage conditions

T = 10-30°C, RH 30-80 %, keep out of direct light.

Transportation conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

5 years.

Packing

Aluminum foil bags for 200 or 500 units.

Packing's information: product code and description, process for intended use, instructions for use, product presentation, classification according to regulation, batch number, manufacturer information and expiration date.

Labelling

On product's label: product code and description, process for intended use, classification according to regulation, batch number, expiration date, bar code and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information

Read product's instructions for use thoroughly before use.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat or any sterilization processes other than Steam.

Do not reuse the sterilizer until the indicator has changed into the color indicated in the product.